

# TECHNICAL DATA SHEET




## Product Name and Code:

Chemical Indicator Label for Steam– **LSL10**

## Product Description:





LISTER Chemical Indicator Label for Steam LSL10 consists of a non-cellulosic plastic backing material with a pressure-sensitive adhesive on one side and chemical indicator spot on the other side. The chemical indicator spot turn from yellow to black after exposure to steam.



## Indications for use:

Use the LISTER Chemical Indicator Label LSL10 as an external pack process indicator to differentiate unprocessed items from items processed in sterilizers. The chemical indicator stripes turn from yellow to black after exposure to steam.

## Performance characteristics:

Reference Color Change Conditions	Illustration
Saturated steam, 121 °C, $\geq 2$ mins, $\leq 10$ mins; Saturated steam, 134 °C, $\geq 0.3$ mins, $\leq 2$ mins;	 UNPROCESSED
	 FAILURES
	
	 PROCESSED AND SATISFACTORY

## Technical information:

Weight	440-450g / bag
Period of validity	36 months
Dimensions	18cm*24.1cm

## Compliance:

ISO 11140-1:2014 Sterilization of health care products-Chemical indicators-Part 1: General requirements.

Manufactured by: LISTER BIOMEDICAL CO.,LTD.  
Building 2, No.72 Shengli Road, Yixiu, Anqing, Anhui, China.  
Tel. +86 556 5689070  
E-mail: [info@listerbiomed.com](mailto:info@listerbiomed.com) -[www.listerbiomed.com](http://www.listerbiomed.com)

***Package:***

13500PCS/BAG

**Version: 20250120**